

**National Institute for Occupational Safety and Health
Evaluation and Testing Branch
Test Data Sheet**



Task Number: TN-26024
Test: Exhalation Resistance Test
Manufacturer: Sleepnet Corporation

Reference No.: CFR 84.172
STP No.: 3

#48012 - Envo Pro Mask Body
#39825 - N99 Filter
#39811 - Envo Pro Filter Cage
#55041 - Holo Headgear

Filter Type: Filter Only

Item Tested: 48012, 39825, 39811, 55041

Sample	Maximum Allowable Resistance (mm of H ₂ O)	Actual Resistance (mm of H ₂ O)	Result
	Exhalation	Exhalation	
1	25	10.2	PASS
2	25	10.2	PASS
3	25	9.9	PASS

Overall Result: PASS

Comments:

Was all equipment verified to be in calibration throughout all testing? Yes No

Signature:

Date: 9/28/2022

Engineering Technician

Form revision: 1.1 18 February 2021

National Institute for Occupational Safety and Health
Evaluation and Testing Branch
Test Data Sheet



Task Number: TN-26024
Test: Inhalation Resistance Test
Manufacturer: Sleepnet Corporation
Item Tested: 48012, 39825, 39811, 55041

Reference No.: CFR 84.172
STP No.: 7

Filter Type: Filter Only

Sample	Maximum Allowable Resistance (mm of H ₂ O)	Actual Resistance (mm of H ₂ O)	Result
	Inhalation	Inhalation	
1	35	11.2	PASS
2	35	11.7	PASS
3	35	11.4	PASS

Overall Result: PASS

Comments:

Was all equipment verified to be in calibration throughout all testing? Yes No

Signature:

Date: 9/28/2022

Engineering Technician

Form revision: 1.1 18 February 2021

Task Number: TN-26024

Reference No.: CFR 84.173

Test: Exhalation Valve Leakage Test

STP No.: 4

Manufacturer: Sleepnet Corporation

Item Tested: 39825


Comments:

Was all equipment verified to be in calibration throughout all testing?

Yes

No

Signature:



Date: 9/29/2022

Engineering Technician